

THE FOCAL SERIES

ALL CARTS INCLUDE:

- Best lock with removable core
- Spring loaded Narcotic box
- Trash compartment
- Pullout writing surface included, additional item on the F297UD model
- Card row dividers with sub-dividers
- Heavy duty caster with brake
- Self closing slides

- Pre-placed inserts for flexible placement of accessories
- Easy change handles for quick color change or refurbishing
- Aluminum shell with vinyl coating insures the ultimate balance of durability and lightweight
- Extended base for maximum stability



Optional double row narcotic drawer with Best® Lock.



F297UD MODEL CAPACITY UP TO 297 UNIT DOSE BOXES

CART SPECIFICATIONS

Cart Base 23" d X 43 5/8" w
 Cart Height 43 1/2"
 Work Surface 21" d X 39" w

DRAWER SPECIFICATIONS

Left Side
 (1) Small Drawer 18" d X 14 1/8" w X 3 1/2" h ID
 (3) Large Drawers 18" d X 14 1/8" w X 9 7/8" h ID
 Right Side
 (9) Small Med Drawers 18" d X 23" w X 3 1/2" h ID

F288UD MODEL CAPACITY UP TO 288 UNIT DOSE BOXES

CART SPECIFICATIONS

Cart Base 24 1/2" d X 34 3/4" w
 Cart Height 43 1/2"
 Work Surface 22 1/2" d X 33 3/4" w
 Slide out shelf 16" d X 13 1/2" w

DRAWER SPECIFICATIONS

(7) Small Drawers
 19 3/8" d X 28 3/8" w X 3 1/4" h ID
 (1) Large Drawer
 18" d X 28 3/4" w X 9 7/8" h ID



F198UD MODEL CAPACITY UP TO 198 UNIT DOSE BOXES

CART SPECIFICATIONS

Cart Base 24 1/2" d X 28" w
 Cart Height 43 1/2"
 Work Surface 22 1/2" d X 27" w
 Slide out shelf 16" d X 13 1/2" w

DRAWER SPECIFICATIONS

(7) Small Drawers
 19 3/8" d X 22 1/8" w X 3 1/4" h ID
 (1) Large Drawer
 18" d X 22 1/8" w X 9 7/8" h ID

OPTIONAL FEATURES


Customize your cart by specifying the following options:



All Purpose Organizer Tray
 Holds cups, straws, etc. Attaches to the side of the cart.
 Place on work surface for easy access.
 Fits perfectly on top of Exterior Tilt Bin Organizer.
Item #01501



**All Purpose Organizer Tray
 with Exterior Tilt Bin Organizer**
 Holds spoons, cups, straws, etc.
Tilt Bin Item #T2001



**New Stor-Flex
 Externals Drawer Divider System**
SMALL F3, APS-36, F198UD Item #E403
MEDIUM F4, APS-48, F288UD Item #E404
LARGE F5, F297UD Item #E405



Silent Knight Pill Crusher
Item #P1611



One and Three Row Exchange Bins
One Row Item #B203 Three Row Item #B204



Sage Sharps Container
 Includes locking security feature.
Item #S1802

Your Choice Of Colored Drawer Handles



Post Office Box 888 El Segundo, California 90245 (800) 437-6633

www.mmimedcarts.com